

Suction Irrigation Set

✓ Control ✓ Comfort ✓ Reliability ✓ Versatility




Trumpet valve

- Trumpet valve design enables fine control of suction and irrigation flows.
- Spring loaded mechanism allows fingertip control and provides good haptic feel with immediate return and shut off.



Comfortable, ergonomic handle

- Comfortable, ergonomic handle design providing high levels of operator precision and control.



Flexible tubing

- Soft, flexible tubing is highly pliable with low memory meaning that it is easily handled with minimal resistance to the operator.
- Anti-kink design and wide lumen of tubing ensure constant and continued performance.

Suction Irrigation Set




Spike configurations

- Single and dual spike configurations
- Laser channel

Removable, interchangeable probes

- 5mm x 33cm
- 5mm x 45cm
- 10mm x 33cm



Suction Irrigation Set	Description	Units per box
LM201	5mm x 33cm Suction Irrigation Set (Single Spike)	Box – 20
LM202	5mm x 33cm Suction Irrigation Set (Dual Spike)	Box – 20
LM2025	5mm x 45cm Disposable Suction Probe	Box – 10
LM20210	10mm x 33cm Disposable Suction Probe	Box – 10



T: +41 791008855
E: Info@perfomed.ch
Place Chauderon 18
1003 Lausanne
Suisse

<https://perfomed.ch/>